

ZipThaw 202™ & ZipSleeve

# ZipFast. ZipSafe. ZipThaw.

ZipThaw is a plasma thawing medical device that is fast, dry, portable, precise, easy-to-use. Most critically, ZipThawing plasma recovers the maximum levels of Coagulation Factors and Antibodies.



## Melting & thawing at the point of care.

### Safe

Disposable ZipSleeve triple redundancy fail-safes dramatically reduces cross-contamination.

### Easy

Integrated barcode scanner, touch screen, disposable sleeve—and a dry & portable form factor make ZipThaw easy to use.

### Preserve Levels

ZipThawing plasma is proven to preserve the maximum levels of coagulation factors and antibodies.

### Fast

Melt or Thaw plasma from frozen in minutes, all at the point of care.

### Portable

Dry, light weight and never needs recalibration when moved.

### Quality

Made in the USA. Fully certified. **3-year full replacement warranty.**

ZipThaw is a dry thawing medical device that lets you quickly, evenly and precisely melt and thaw fresh frozen plasma.

Use the integrated barcode scanner, slide a plasma specimen bag in the disposable ZipSleeve, insert into ZipThaw chamber, press a button for Melt or Thaw presets, stop the cycle anytime.

- ▶ **Accelerate throughput:** Each chamber works independently so you can continuously melt and thaw specimens.
- ▶ **Know viability:** Precise temperature display with alarm that alerts shut off.
- ▶ **Capture data:** Each thaw is tracked on a hard drive with data to export later.

**Contact us to learn more about ZipThaw.**

Web: [fremonscientific.com](http://fremonscientific.com)

Email: [info@fremonscientific.com](mailto:info@fremonscientific.com)

Phone: 1-617-308-4405

## ZipThaw 202™ & ZipSleeve

# Maximum Preservation of Coagulation Factors and Antibodies

- ▶ In 2019, at the San Diego Blood Bank and UC San Diego School of Medicine, ZipThaw reproducibly and reliably thawed FFP and PF24 with preservation of coagulation factors for clinical use - data submitted to FDA and presented at ASCO and AABB 2019 conferences.
- ▶ In 2020, ZipThaw demonstrated to thaw COVID-19 convalescent plasma to reliably recover immunoglobulins: IgA, IgM, and IgG mg/dl levels were measured across all samples by the UCSD Clinical Chemistry Lab. 100% concordance with all samples yielding equivalent levels of Immunoglobulins recovered post-ZipThaw compared to pre-freeze levels.



ZipSleeve is a disposable protective layer with patented sensors. It works with ZipThaw to continuously monitor and report your specimen's temperature, not its surroundings. It is a multi-use disposable that tracks its use status and can't be reused past expiration.

## Specifications

**Thawing Technology:** Electrical power, disposable RFID sleeve

**Typical Capacity:** Single or dual independently operated chambers

**Heating Temperature:**

- Thawed plasma not to exceed 37°C±0.5°C
- Precise sensors display end-of-run actual thawed plasma temperature

**ZipSleeve Sensors Cutoff Temperature:** Melt 15°C / Thaw 31°C

**Agitation Method:** Electronic control mechanical massage

**Visual Display:** 7" touch screen

**System:** Main CPU: VAR-SOM-SOLO/Dual: Freescale i.MX6

**Storage humidity (non-condensing):** 30% - 90% non-condensing

**Interface & Network:** USB 2.0 OTG

**Internal Memory:** Stores 2,000 most recent cycles

**Internal Power Supplies:** Medical grade, isolating power-supply

- Input: 100-240Vac, maximum rated current 2.5A
- Outputs: 2 X 24VDC, 12.5A maximum current

**Electrical Rating (V/Hz/A):** 100-240Vac, 50/60Hz 2.5Amp

**Size: [W x D x H]:**

- cm: 47.1 x 27.6 x 41.6 (including the top handle)
- inch: 18.5 x 10.9 x 16.4 (including the top handle)

**Weight:** 29lb/13Kg

**Portable:** Yes

**Environmental Requirements:**

- Ambient Operating Temp: 10°C to 32°C (50°F to 86°F)
- Ambient Humidity: 20% to 70% non-condensing
- Storage Temperature (in original packaging): -20°C to 70°C (-4°F to 158°F)